



## ECOJECT® Plus Safety syringes & NEOJECT® Plus Safety needles

Dispomed® safety system.

ECOJECT® SAFETY disposable syringe in connection with NEOJECT® PLUS SAFETY offers an effective protection from needlestick injuries and the risk of infection in accordance with TRBA250.

- 2-part syringe Luer-Lock and needle, for single use
- Ideal readability due to high contrast black graduation
- Peel open package
- Sterile EO (ethylene oxide) in accordance with ISO 11135-1
- Product classification IIa

[Here you can download the instruction of use](#)



### Ecject® Safety 2-part disposable syringe

	Art.-No.	Pack. unit
2-3 ml	20202	100 / 3000
5-6 ml	20205	100 / 2500
10-12 ml	20210	100 / 1500

### Neoject® Plus safety needles

	Art.-No.	Pack. unit
No.1 (0,90 x 40mm) 20G x 1 1/2"	11001	100 / 3000
No.2 (0,80 x 40mm) 21G x 1 1/2"	11002	100 / 3000
No.12 (0,70 x 30mm) 22G x 1 1/4"	11012	100 / 3000
No.14 (0,60 x 30mm) 23G x 1 1/4"	11014	100 / 3000
No.16 (0,60 x 25mm) 23G x 1"	11016	100 / 3000