



## ANEL Special hypodermic needles

Special blunt needles for the optical and dental purpose

- blunt needles for single use
- straight or curved
- sterile R (Gamma)
- precision lancet grinding
- dry-siliconization
- transparent Luer / Luer-Lock tip
- protection cap
- small dead space
- exact glueing



				<b>Art.No.</b>	<b>pack. unit</b>
straight	0,40 x 18	27G x 3/4"	grey	10090	100 / 3000
straight	0,50 x 18	25G x 3/4"	orange	10091	100 / 3000
straight	0,60 x 18	23G x 3/4"	blue	10092	100 / 3000

### curved in 14 mm distance from tip (45° angle)

				<b>Art.No</b>	<b>pack. unit</b>
curved	0,40 x 18	27G x 3/4"	grey	10093	100 / 3000
curved	0,50 x 18	25G x 3/4"	orange	10094	100 / 3000
curved	0,60 x 18	23G x 3/4"	blue	10095	100 / 3000