


**Dispomed Witt oHG**  
Am Spielacker 10 -12, D-63571 Gelnhausen  
Tel.: (+49) 0 60 51 - 96 23 -0  
Fax.: (+49) 0 60 51 - 96 23 -0  
eMail: info@dispomed.de



**Dispomed®**

## DISPOMED® bottle-needles

- bottle-needle for single use
- non-coring bevel
- dry siliconisation
- high flow rate
- coloured hub, ISO colour code
- sterilizable by gamma radiation (max. 25kGy)
- sterilizable by Ethylenoxide (max. 60°C)



model code	colour	cannula	Pack. unit
13010	white	Ø1,04 x 30mm	1000 / 5000
13007	pink	Ø1,22 x 30mm	1000 / 5000
13011	pink	Ø1,22 x 40mm	1000 / 5000
13008	blue	Ø1,53 x 25mm	1000 / 5000
13005	red	Ø1,83 x 35mm	1000 / 5000
13006	green	Ø2,07 x 35mm	1000 / 5000