



Dispomed® disposable syringe Insulin

3 part syringes without dead space

- according to DIN EN ISO 8537
- LUER, centric
- no loss of pharmaceuticals by displaced piston
- smooth running piston
- red scale for a better reading
- sterile EO (Ethylene oxide)



			Art.-No.	pack. unit
1ml Insulin U-40	without	needle	22001	100 / 2000
1ml Insulin U-100	without	needle	22002	100 / 2000
1ml Insulin U-40	0,45 x 13	26G x 1/2"	22003	100 / 2000
1ml Insulin U-40	0,40 x 13	27G x 1/2"	22033	100 / 2000
1ml Insulin U-40	0,33 x 13	29G x 1/2"	22023	100 / 2000
1ml Insulin U-100	0,45 x 13	26G x 1/2"	22004	100 / 2000
1ml Insulin U-100	0,40 x 13	27G x 1/2"	22005	100 / 2000
1ml Insulin U-100	0,33 x 13	29G x 1/2"	22043	100 / 2000