



Dispomed® Disposable syringes TBC / Heparin

3 part syringes without dead space

- LUER, centric
- Latex- and PVC-free
- no loss of pharmaceuticals
- smooth-running, acytotoxic piston
- for TBC: black scale, indelible graduation: 0,01 ml
- for Heparin: red scale, indelible graduation: 0,05 ml
- sterile EO (Ethylene oxide)



			Art.No.	pack. unit
1 ml TBC	without	needle	22006	100 / 2000
1 ml TBC with needle	0,50 x 16	25G x 5/8"	22007	100 / 2000
1 ml Heparin 25000 I.U.	without	needle	22009	100 / 2000
1 ml Heparin 25000 I.U. with needle	0,50 x 16	25G x 5/8"	22010	100 / 2000