



## ECOJECT® PLUS

Economical disposable syringes according to DIN EN 7886

- **Latex- and PVC-free**
- sterile EO (Ethylene oxide)



### ECOJECT® PLUS Perfusion

Pump syringe 50ml with mark on 60ml fit for Braun Perfusor and compatible devices

	Art.No.	pack. unit
Luer-Lock, transparent	21051	50 / 300
Luer-Lock, transparent, with needle 1,8 x 25 mm	21052	50 / 300
Luer-Lock, <b>orange coloured</b> for lightsensitive pharmaceuticals with attached needle 1,8 x 25 mm	21064	50 / 300

### ECOJECT® PLUS Infusion

Pump syringe 50ml with mark on 60ml fit for Braun Perfusor and compatible devices

	Art.-No.	pack. unit
Luer-Lock, transparent	21050	50 / 300
Luer-Lock, transparent, with needle 1,8 x 25 mm	21054	50 / 300

### ECOJECT® PLUS Luer-Syringe 50 ml, with 60 ml graduation

- 50/60 ml syringe, 3part
- LUER, excentric
- Graduation: 2ml
- Latexfree
- EO sterilized



	Art.-No.	pack. unit
50/60 ml Luer-Syringe	21500	50 / 50