



## INFUJECT® Pump syringes

30 ml and 50/60 ml syringes, 3 part, adapted for Fresenius Injectomat and compatible devices

- According to DIN EN ISO 7886
- Luer Lock, centric
- Pressure tested to 6 bar
- Steril R (Gamma sterilized)
- Rub firm graduation
- Cylinder and plunger made of Polypropylene
- latexfree
- pistonstop
- exact reading of the syringe- volume by clearly 1 ml graduation
- for the application of photosensitive pharmaceuticals: with orange (amber) colouring
- sterile R (Gamma radiation)



				<b>Art.-No.</b>	<b>pack. unit</b>
30 ml with blanking plug, <b>without needle</b> , transparent, "Heparinsyringe"				22030	100 / 600
50/60 ml with blanking plug, <b>without needle</b> , transparent,				22050	50 / 300
50/60 ml, transparent, <b>with</b> attached needle	2,0 x 30mm	14G x 1 1/4"		22054	50 / 300
50/60 ml, orangefarben, <b>with</b> attached needle	1,8 x 25mm	15G x 1"		22064	50 / 300