



NEOJECT® Hypodermic needles, special sizes

Extremely thin walled needles for the intravenous, intramuscular, subcutaneous oder intraarteriel injection

- Precision lancet grinding, long bevel
- Dry siliconisation
- Exact glueing
- High flow rate
- Transparent luer-lock hub, ISO colour code 6009
- Peel-open package
- Sterile EO (ethylene oxide) in accordance with ISO 11135-1



				Art.-No.	Pack. unit
Insulin	0,33 x 13	29G x 1/2"	red	10083	100 / 3000
Insulin	0,42 x 13	27G x 1/2"	grey	10127	100 / 3000
Insulin	0,45 x 13	26G x 1/2"	brown	10027	100 / 3000
TBC	0,50 x 16	25G x 5/8"	orange	10026	100 / 3000
Subcutaneous	0,50 x 25	25G x 1"	orange	10089	100 / 3000
Intramuscular	0,60 x 60	23G x 2 3/8"	blue	10025	100 / 3000
Intramuscular	0,70 x 25	22G x 1"	black	10036	100 / 3000
Intramuscular	0,70 x 40	22G x 1/2"	black	10021	100 / 3000
Intramuscular	0,70 x 50	22G x 2"	black	10022	100 / 3000
Intramuscular	0,80 x 25	21G x 1"	green	10024	100 / 3000
Intramuscular	0,80 x 50	21G x 2"	green	10028	100 / 3000
Intramuscular	0,90 x 25	20G x 1"	yellow	10023	100 / 3000
Intramuscular	0,90 x 50	20G x 2"	yellow	10029	100 / 3000
deep intramuscular	0,90 x 70	20G x 2 3/4"	yellow	10030	100 / 3000
Blood withdrawal	1,10 x 25	19G x 1"	ivory	10053	100 / 3000
Blood withdrawal	1,10 x 30	19G x 1 1/4"	ivory	10031	100 / 3000
Blood withdrawal	1,10 x 40	19G x 1 1/2"	ivory	10032	100 / 3000
Blood withdrawal	1,10 x 50	19G x 2"	ivory	10070	100 / 3000
Blood withdrawal	1,20 x 40	18G x 1 1/2"	pink	10033	100 / 3000
Blood withdrawal	1,20 x 50	18G x 2"	pink	10068	100 / 3000
Veterinary	1,50 x 50	17G x 2"	ivory	10035	100 / 3000
Veterinary	1,60 x 25	16G x 1"	white	10077	100 / 3000
Veterinary	1,60 x 50	16G x 2"	white	10076	100 / 3000
Sample needle	1,80 x 40	15G x 1 1/2"	blue-grey	10034	100 / 3000