



NEODENT® dental hypodermic needle

Disposable needles according to
DIN EN ISO 7864

- precision lancet grinding
- dry-siliconization
- transparent Luer / Luer-Lock tip, protection cap
- Colour code according to DIN EN ISO 6009
- small dead space
- exact glueing
- high flow rate
- sterile R (Gamma radiation)



				Art. no.	Pack. unit
dental injection	0,50 x 23	25G x 7/8"	orange	12001	100/3000
dental injection	0,52 x 42	25G x 1 5/8"	orange	12002	100/3000
dental injection	0,42 x 23	27G x 7/8"	grey	12003	100/3000
dental injection	0,42 x 23	27G x 1 5/8"	grey	12004	100/3000
blunt dental injection	0,35 x 25	23G x 1"	blue-green	12005	100/3000