



DISPOMED® Injectstoppers

- safety closing stoppers
- membrane Latex free!
- with injection port for pharmaceuticals
- Luer-Lock connector, female
- sterile EO (Ethylene oxide)



	Art. no.	Pack. unit
Injectstoppers, yellow	70013	100 / 200

DISPOMED® Universal stoppers

- safety closing, no injection port
- double function:
Luer Lock male and Luer Lock female
- pressure resistant up to 200kPa (2 bar), universal use for pressure and gravity systems
- Steril EO (Ethylene oxide)



	Art. no.	Pack. unit
Kombistopfen, orange	70014	100 / 200
Kombistopfen, red	70015	100 / 200
Kombistopfen, blue	70016	100 / 200