

**Dispomed Witt oHG**  
Am Spielacker 10 -12, D-63571 Gelnhausen  
Tel.: (+49) 0 60 51 - 96 23 -0  
Fax.: (+49) 0 60 51 - 96 23 -0  
eMail: info@dispomed.de



**Dispomed®**

## DISPOMED®Pleura Puncture Set

### Pleura Puncture Set

completely closed system composed of:

- 50 ml syringe(Luer Lock) for single use
- 3 way stop cock
- tube 30 cm , duplex LuerLock - male
- pleura puncture needle 1,8 x 80mm
- secretion bag 1500ml or 2000ml
- safety package (double blister)
- sterile EO (Ethylene oxide)



	<b>Art.-No.</b>	<b>pack. unit</b>
DISPOMED® Pleura-Puncture Set, 1500 ml	80050	40 / 40
DISPOMED® Pleura-Puncture Set, 2000 ml	80051	40 / 40