



HEMOMED® Transfusionsets

Blood administration set for single use acc. to IN 58360-1 und ISO 1135-4

- Long drip chamber
- 200 µm bloodfilter
- Injection port: latexfree
- Precision drip controller
- Luer and Luer-Lock fitting
- Standard tube 150cm
- Spike for bag and bottle
- Sterile R (Gamma radiation)



		Art.-No.	pack. unit
With ventilation valve	BLL	60003	100 / 100
Without ventilation valve	LL	60011	100 / 100
Without ventilation valve; without injection port	LLP	60111	100 / 100
#60111	Suitable for use with blood bag pumps		

P: pressure tested up to 200 kPa (2 bar).

suitable for use with blood bag pumps (color code of drip control: orange)

L = Luer LL = Luer Lock P = pressure B = ventilation valve N = needle