



PERFULINE® Infusion tubes for syringes pumps, PVC-free

PERFULINE® infusion tubes acc. to ISO 8536 / DIN 58362-SL-P

- Ecological , because of PVC-free material
- pharmaceuticals are completely available, no ad- or absorption
- Luer-Lock and Luer fitting; with protection cap
- compatible for all equipments with authorization of standard infusion tubes acc. to DIN/ISO (e.g. Braun Perfusor)
- leakage tested, pressure of 500 kPa (5 bar)



material				Art.-No.	pack. unit
PE	1,0 x 2,0	150cm	transparent	50155	100 / 200
PVC-free	1,5 x 2,6	150cm	transparent	50160	100 / 200
PVC-free	1,5 x 2,7	200cm	transparent	50202	100 / 200
for light sensitive pharmaceuticals:					
PE	1,0 x 2,0	150cm	amber	50166	100 / 200