



SOLUFLO® Infusionsets for operation under pressure (P) and gravity (G)

Infusion set for single use
 in accordance with ISO 8536-8 P (Pressure)

- Spike for bag and bottle
- Air filter
- Fluid filter: 15µm pore size; 120mm² filter area
- Long drip chamber
- Tubing: outside diameter 4,10mm; inside diameter 3,00mm
- Precision roller clamp
- All product variations are latex free
- Sterile R (Gamma radiation) in accordance with ISO 11137-1



Product variations:

L = Luer

LL = Luer Lock

P = for operation under pressure (colour of roller clamp: orange)

V = for non-return valve (colour of the rollerinside of the roller clamp: orange; Protective cap for Luer Lock adapter: orange)

Item No.	type	tube length	inject. port		pack. unit
60012	LLP	150cm		PVC (DEHP-free)	100 / 100
60913	LLP	190cm		PVC (DEHP-free)	100 / 100
60014	LLVP	150cm		PVC-free	100 / 100
61214	LLVP	175cm		PVC (DEHP-free)	100 / 100
60214	LLVP	180cm		PVC (DEHP-free)	100 / 100
60006	YLLP	180cm	Y-port	PVC (DEHP-free)	100 / 100
60618	YLLP	180cm	Y-port	PVC (DEHP-free)	100 / 100