



SOLUFLO® Infusion set for operation under Gravity (G)

Infusion set for single use in accordance with ISO 8536-4 G (Gravity)

- Spike for bag and bottle
- Air filter
- Fluid filter: 15µm pore size; 120 mm² filter area
- Long drip chamber
- Tubing: outside diameter 4,10mm; inside diameter 3,00 mm
- Precision roller clamp (blue)
- All products variation are latex free
- Sterile R (gamma radiation) in accordance with ISO 11137-1

Product variations

L = Luer

LL = Luer Lock

Vet = for veterinary



Item No.	Type			Tubing	Pack.unit
60009	L	with injection port	PVC-free	150 cm	100 / 100
60005	LL	with injection port	PVC-free	150 cm	100 / 100
60518	LL	with injection port	PVC-free	180 cm	100 / 100
60018	LVet	drip chamber without fluid filter; with injection port	PVC-free	150 cm	100 / 100
60017	LLVet	drip chamber without fluid filter; with injection port	PVC-free	150 cm	100 / 100
60055	LL	with injection port	PVC (DEHP-free)	150 cm	100 / 100
60619	LL	Soluflo: Item 60005 with Ecoflo Safety Scalp Vein Set 19G x 3/4" (1,1x19mm)	60005: PVC-free Ecoflo: PVC (DEHP-free)	150 cm; 10 cm	100 / 100