



Dispomed® Straussneedles

- Ultra thin wall needle acc. to DIN EN ISO 9626
- hub and gripplate: Polypropylen
- gripplate: checkered and non slipped
- Protection hub: Polyethylen
- Precision lancet bevel
- Dry siliconisation
- Clean bonding
- Peel open package
- Sterile R (Gamma radiation)
- Luer



| | | | Art.-No. | pack. unit |
|-----------------|--------------|--------|-----------------|-------------------|
| needle 0,9 x 43 | 20G x 1 3/4" | yellow | 13009 | 100 / 1000 |
| needle 1,2 x 43 | 18G x 1 3/4" | pink | 13012 | 100 / 1000 |
| needle 1,5 x 43 | 17G x 1 3/4" | blue | 13015 | 100 / 1000 |
| needle 1,8 x 43 | 15G x 1 3/4" | orange | 13018 | 100 / 1000 |
| needle 2,0 x 4 | 14G x 1 3/4" | green | 13020 | 100 / 1000 |

unsteril / bulk on request