



VASUFLO®-Safety scalp vein set

- scalp vein set with active protection clip
- after use: move the safety protector over the needle
- minimizing the injury risk for user and patient
- the scalp vein sets are available in the established gauge sizes with 30cm tubing
- needle: precision grinding, dry siliconisation of the inner and outer side, inclination 8°
- safety tubing: highly elastic, buckling resistant
- coloured size marking
- luer lock fitting with cap



GAUGE	Ø	length	item no.	pack. unit
18 G	1,20 mm	19 mm	46018	100 / 1000
19 G	1,10 mm	19 mm	46019	100 / 1000
20 G	0,90 mm	19 mm	46020	100 / 1000
21 G	0,80 mm	19 mm	46021	100 / 1000
22 G	0,70 mm	19 mm	46022	100 / 1000
23 G	0,60 mm	19 mm	46023	100 / 1000
25 G	0,50 mm	19 mm	46025	100 / 1000
27 G	0,40 mm	19 mm	46027	100 / 1000

unsteril / bulk or other tubing lengths on request

for example: 10cm tubing, sterile