



## VASUFLO®INT intravenous catheters

Peripheral i.v.-catheter with injection port acc. to DIN EN ISO 10555-5

- thin walled Teflon (FEP) catheter for high flow rate
- excellent fixing on the skin
- colour code acc. ISO standart
- longitudinal injection port
- Luer-Lock fitting
- special bevel cut for minimal puncture trauma
- Precicely fitting and colour coded Mandrins



|                  |              |     |        | <b>Art.-No.</b> | <b>pack. unit</b> |
|------------------|--------------|-----|--------|-----------------|-------------------|
| needle 2,0 x 45  | 14G x 1 3/4" | 14G | orange | 42014           | 50 / 500          |
| needle 1,75 x 45 | 16G x 1 3/4" | 16G | grey   | 42016           | 50 / 500          |
| needle 1,5 x 45  | 17G x 1 3/4" | 17G | white  | 42017           | 50 / 500          |
| needle 1,3 x 45  | 18G x 1 3/4" | 18G | green  | 42018           | 50 / 500          |
| needle 1,3 x 32  | 18G x 1 1/4" | 18G | green  | 42019           | 50 / 500          |
| needle 1,0 x 32  | 20G x 1 1/4" | 20G | pink   | 42020           | 50 / 500          |
| needle 0,9 x 25  | 22G x 1"     | 22G | blue   | 42022           | 50 / 500          |
| needle 0,74 x 19 | 24G x 3/4"   | 24G | yellow | 42024           | 50 / 500          |
| needle 0,64 x 19 | 26G x        | 26G | violet | 42026           | 50 / 500          |