




## VASUFLO® Scalp vein set

Device for vein puncture in periphere parts and connection for injection and infusion therapy, short term use (<5 days)

- Needle: precision lancet bevel; siliconisation of inner and outer side; inclination 8°
- Safety tubing; non buckling; elastic; 30cm
- Colour-coded butterfly
- Luer Lock fitting with cap
- Sterile R (Gamma radiation)



			<b>Art.-No.</b>	<b>pack. unit</b>
needle 1,2 x 19	18G x 3/4"	pink	40018	100 / 1000
needle 1,1 x 19	19G x 3/4"	ivory	40019	100 / 1000
needle 0,9 x 19	20G x 3/4"	yellow	40020	100 / 1000
needle 0,8 x 19	21G x 3/4"	green	40021	100 / 1000
needle 0,7 x 19	22G x 3/4"	black	40022	100 / 1000
needle 0,6 x 19	23G x 3/4"	blue	40023	100 / 1000
needle 0,5 x 19	25G x 3/4"	orange	40025	100 / 1000
needle 0,4 x 19	27G x 3/4"	grey	40027	100 / 1000

unsterile / bulk or other tubing lengths: on request