




DISPOMED® Biopsy Needles

- for soft tissue biopsy
- for biopsy of tissue samples
- 20 mm recessed specimen notch
- Triple bevel "V" cannula - clean cut
- Stainless steel, high rigidity
- grip comfortable fits to finger, the ergonomic design ensures a precise puncture
- Easy handling because of the "lock"-automatic
- sterile EO (Ethylene oxide)



		Art. no.	pack. unit
2,10 x 75 mm	14G x 3"	70007	1 / 50
2,10 x 114 mm	14G x 4 1/2"	70008	1 / 50
2,10 x 152 mm	14G x 6"	70009	1 / 50
2,10 x 200 mm	14G x 8"	70010	1 / 50
1,60 x 75 mm	16G x 3"	70022	1 / 50
1,60 x 114 mm	16G x 4 1/2"	70023	1 / 50
1,60 x 152 mm	16G x 6"	70024	1 / 50
1,20 x 75 mm	18G x 3"	70025	1 / 50
1,20 x 114 mm	18G x 4 1/2"	70026	1 / 50
1,20 x 152 mm	18G x 6"	70027	1 / 50